



"KOHARU" - Elderly Care Manikin

● Accessories/Pajamas, Baby powder
<Wig sold as a separate item>

For elderly care training
in good relations



Wig (extra charge)



Reproduction of the bending of the back peculiar to the elderly.

In order to make it similar to the movements of the elderly, the back is made so that it can bend forward, making realistic training possible such as when seating in a wheelchair.

The movement of legs and hip joints are smooth.

Forward/backward and outward movement of the legs are more practical, and the hip joints open widely. Realistic training in helping excretion using an excrement or toilet is possible.

Joints are made subtle, hardly noticeable.

Since it is completely seamless, training in bathing can be done safely and smoothly.

Durable yet easy-to-handle skin material.

The skin is made of high-quality silicone, and can be used freely with a shower, sponge, or detergent.

Feature

1. The head moves in all directions (front, back, right, left).



The head moves not only front and back, but also right and left, making practical training possible as if on a real human body.

2. Movement of the joints are more real.



A bone structure is built-in. Each joint in the shoulders, elbows, wrists, hips, thighs, knees, and ankles move rather stiffly as in the care of actual elderly people. Training in exchanging pajamas and bedding, bathing, and cleaning can be done more realistically.

Specifications

Height	150cm
Weight	About 12kg
Material	Silicone

Practice

Oral hygiene care



The denture is removable, and care of the oral cavity can be practiced.

Excretion and urination



It is possible to practice not only assistance in excretion using an excrement or toilet, but also in removing excrement from the anus.

Wheelchair assistance



The joints in the hips, knees, and ankles bend, making it possible to learn smoothly the way to offer wheelchair assistance.

Changing position



Due to the new function in the joints of each body part, assistance in getting up to prevent bedsores, balance in sitting, and changing position can be offered more realistically.

Bathing



When practicing bathing (whole or partial), it is possible to put in actual hot water. Can also use a shower, sponge, or detergent.