


M100

"TOMOKO" - Nursing Training Manikin

● Accessories/Wig for shampoo, Pajamas, Catheters

Mate for nursing training

- Average size of adult woman
- Reasonable weight for practice
- Can feel the softness of skin like treating real patient
- Lid of abdomen and chest fit perfectly
- Continuous bladder catheterization is available with reality (indwelling/ balloon catheterization can be performed)
- Condition of catheter can be checked from out side



Practice



- Stoma care
- Practice with balloon catheter available
- Bed bath
- Shampoo
- Posture change
- Transportation
- Changing pajamas
- Nutrition (nasal/oral)
- Assisting excretion
- Restriction
- Stoma care
- Injection
- Bandaging
- Fomentation
- Enema
- Manual disimpaction
- Catheterization
- Gastric lavage
- Intestine lavage
- Tracheotomy care

Specifications

Height	158cm
Weight	About 15kg
Material	Soft vinylchloride

M100-1 Interchangeable to male genitalia model

Male Genitalia Model for "TOMOKO" (M100)

● Weight/1 piece 0.5kg ● Accessories/Irrigator, Intestines



*Irrigator stand is not included.


M100-2

"KIYOKO" - Bathing Care Manikin

● Accessories/Pajamas
<Wig and Male genitalia model sold as separate items>

Various joints move naturally as those of human body, and can simulate with feeling close to actual patient.

- The lower back smoothly bends anteriorly and posteriorly. The model can be seated in a wheelchair for simulation with a feeling similar to seating a real patient.
- The head can be turned laterally to 30° in both directions as well as upward and downward so that the model can take nearly the same postures as the human body.



- Since the lower limbs can be freely moved back and forward or abducted, the legs can be spread sufficiently. Practice of assistance in excretion using a bedpan or a urinal can be performed with a true to life feeling.
- Various joints including the shoulders, elbows, wrists, hips, knees, and ankles move naturally as those of the human body and enable smooth execution of practice of changing clothes and bedding.

Practice

- Can actually be bathed.
- Water spontaneously flows into the model so that it can be dipped into water. Drained quickly when it is lifted up from water.
- Tolerates the use of showers, sponge, and detergents.



Specifications

Height	160cm
Weight	About 13kg
Material	Vinyl chloride