

M107-1

"Mamoru"

Bathing model (boy)

● Accessories/Umbilical cord, Baby clothes

M107-2

"Haruko"

Bathing model (girl)

● Accessories/Umbilical cord, Baby clothes

M107-3

SAKAMOTO Baby (boy)

● Accessories/Umbilical cord, Baby powder, Baby clothes

M107-4

SAKAMOTO Baby (girl)

● Accessories/Umbilical cord, Baby powder, Baby clothes

The durability baby model for maternal nursing and maternal classes



Specifications

Height	45cm
Weight	About 3kg
Material	Soft vinyl chloride



Umbilical cord (fit-on & off available)

- The limbs have seamless structures.
- The neck turns sagittally and horizontally. (Please hold the model by supporting the head as you would do in holding a neonate.)
- The weight of the head has the same proportion as that of a neonate.
- The anterior fontanelle and ribs can be palpated.
- Water infused into the mouth flows out from the urethra so that changing diapers can be practiced.

Practice

- Bathing ● Observation measurements ● Umbilical treatment ● Anal thermometry ● Baby exercise ● Changing diapers ● Dressing and undressing



Specifications


Height	49cm
Weight	About 3.3kg
Material	Silicone

Practice

1. Bathing



2. Baby exercise and feeding



3. Dressing / undressing and changing diapers



4. Anal thermometry and various measurements



Others Usable also for the practice of umbilical treatment, cleaning of the acoustic meatuses and nostrils, and guidance at infant screenings.

M125-1

Hospital Training Baby

Practice

- Temporal artery puncture (Temporal artery is exchangeable.)
- Catheter can be put in stomach through transnasal and oral (listen to sounds).
- Catheterization ● Enema ● Injection to the thigh and buttocks. (Injection pads are exchangeable.)
- Stoma ● Umbilical view catheter

Specifications

Height	64cm
Weight	About 3kg
Material	Soft Synthetic resin



4 weeks baby

Boy and girl can be changed by exchanging abdominal cover.

Boy : Abdominal wall with stoma, urethra and bladder.

Girl : Abdominal wall with umbilical view, Urethra and bladder.