

M166-1



Fit-on Female Catheter Simulator

- Weight/About 5kg ● Case size/W55×D26×H27cm
- Accessories/Attachable waterproof pants, Receiving stand (desk-top), Irrigator, Irrigator stand, Catheter, Lubricant oil (Silicone oil), Baby powder, Storage case

Training in catheterization and cleansing of private parts possible.



M166-2



Fit-on Male Catheter Simulator

- Weight/About 5kg ● Case size/W55×D26×H27cm
- Accessories/Attachable waterproof pants, Receiving stand (desk-top), Irrigator, Irrigator stand, Catheter, Lubricant oil (Silicone oil), Baby powder, Storage case



Practice

Practice while wearing



Practice using receiving stand



Training in actual posture

One can wear the attachable pants on a bed, enabling training in actual posture.

Indispensable teaching aid for catheter

It is fitting for practicing as a couple or group.

Realistic learning experience

The external genitals closely resemble the actual object, enabling realistic training in catheter.

Easily attachable on top of clothes

The attachable pants are made of waterproof material and the seams have been waterproof processed, making it possible to wear on top of clothes and practice.

Training in cleansing of private parts possible

As it is waterproof processed, by wearing the attachable pants and setting the genitals, can practice cleansing of private parts.

Can study on a table

By fixing to exclusive "receiving stand", can easily study on a table.

Practice while wearing



Practice using receiving stand

