

M115

Breast Self-Examination Teaching Model

● Case size/W36×D27×H17cm ● Accessories/Baby powder, Storage case




A simulator made of a skin-touch material optimal for learning the self-examination skill.

- Complete self-palpation and observation of the breasts and armpits are possible.
- The true-to-life material makes the model ideal for learning and practicing the self-examination skill of the breasts.
- In one breast, several hardened areas are implemented in regions liable to breast cancer. The other breast is made normal for comparison.

Practice

Practice of palpation for signs of breast cancer

One breast has hardened areas, and the other breast is made normal.



Signs of breast cancer

Lymph nodes.....	1 2
Hardening 1.5mm.....	3 4
Hardening 2.0mm.....	5
Discoloration or depression of the skin.....	6
Pig skin (a state resembling orange peel).....	7
Collapsed nipple.....	8

M114

Model for Breast Massage

● Size/W41×D41×H26cm ● Accessories/Storage bag

- The breast parts, filled with a special resin, give a feeling of real breasts and are appropriate for practicing massage techniques.
- Practice can be performed in the recumbent position due to the wooden stand attached to the back as well as in the upright position.
- With soft nipples, the model is ideal for practicing breast massage shortly after parturition.



M114-1

Mamma Simulator - Breast massage model

● Case size/W48×D48×H18cm ● Accessories/Storage case

The incorporated ribs make the model more real.

- Composed of thoracic ribs and the skin and subcutaneous fat made of special resins feeling like the human body, the model allows practice of puerperal breast treatment and massage with a sense of reality.
- The breast fat in the interior moves according to the hand movement during massage, and the fingers can be thrust into the base of the breast.
- With soft nipples, the model is optimal for practicing breast massage shortly after parturition.

