## M115

# Breast Self-Examination Teaching Model

● Case size/W36×D27×H17cm ● Accessories/Baby powder, Storage case



A simulator made of a skin-touch material optimal for learning the self-examination skill.

- Complete self-palpation and observation of the breasts and armpits are possible.
- The true-to-life material makes the model ideal for learning and practicing the self-examination skill of the breasts.
- In one breast, several hardened areas are implemented in regions liable to breast cancer. The other breast is made normal for comparison.


### Practice

#### Practice of palpation for signs of breast cancer

One breast has hardened areas, and the other breast is made normal.

#### Signs of breast cancer

Lymph nodes······ 1	2
Hardening 1.5mm····· 3	4
Hardening 2.0mm·····	5
Discoloration or depression of the skin·····	6
Pig skin (a state resembling orange peel)······	7
Collapsed nipple	8



# M114-1

## Mamma Simulator-Breast massage model

● Case size/W48×D48×H18cm ● Accessories/Storage case

### The incorporated ribs make the model more real.

- Composed of thoracic ribs and the skin and subcutaneous fat made of special resins feeling like the human body, the model allows practice of puerperal breast treatment and massage with a sense of reality.
- The breast fat in the interior moves according to the hand movement during massage, and the fingers can be thrust into the base of the breast.
- With soft nipples, the model is optimal for practicing breast massage shortly after parturition.

