# Transparent Female Catheter Model

- Accessories/Irrigator, Irrigator stand, Catheter, Lubricant oil, Baby powder, Storage case

## Transparent hypogastrium, inside of bladder observable






Can observe indwelling condition of balloon catheter and flow of urine!



Can practice proper catheterization.

#### Feature

- The external genitalia is made true to nature. You can open the nympha, and the mucous plica etc. are very similar to the real thing.
- It is made so that you cannot find the urethra without opening the nympha.
- The pelvis is built-in. You can practice learning the catheter angle according to the bladder's position.
- You can learn while observing balloon catheter insertion skills and indwelling conditions. You can learn actual resistance and pressure. (Always use a lubricant when inserting catheter.)
- Urine (water) actually flows by catheterization.



# Transparent Male Catheter Model

- Accessories/Irrigator, Irrigator stand, Catheter, Lubricant oil, Baby powder, Storage case



### Practice



Can observe indwelling condition of balloon catheter and flow of urine!



Can practice proper catheterization.

#### Feature

- The external genitalia is made true to nature. You can learn actual resistance and pressure.
- The pelvis is built-in. You can practice learning the catheter angle according to the bladder's position.
- You can learn while observing balloon catheter insertion skills and indwelling conditions. (Always use a lubricant when inserting catheter.)
- Urine (water) actually flows by catheterization.

