


M169

Fit-on Manual Disimpaction Simulator

- Weight/About 4kg ● Case size/W55×D26×H27cm
- Accessories/Belt, Simulated excrement, Lubricant jelly, Baby powder, Storage case
- Director/Hiromi Tashiro

The feel of the anal and rectal areas is realistic.



Practice

Set

Simulated excrement can be easily set and used repeatedly.



Fit-on

By fastening the belt, fitting-on can be easily done.



Wash

The anal area can be removed for cleaning.



Exchangeable parts

- Simulated excrement
- Lubricant jelly