

Postoperative Care Suit

- Weight/About 1.2kg ● Case size/W46×D16×H26cm
- Accessories/Sham tracheal cannula (1 piece), Electrodes (6 pieces), Stoma (1 piece), Suture seal (1 sheet), Male external genitalia (1 piece), Female external genitalia (1 piece), Storage case

A patient is there simply by attaching accessories.

Practicing postoperative care using a naked model is awkward. To cope with such an opinion from the scenes of nursing practice, the postoperative care suit, which can be worn over regular clothes, was developed. The body of the suit, which is waterproof, has holes for abdominal and thoracic drainage, and is packaged with suture seals, a stoma, electrodes, a tracheal cannula, and male and female external genitalia. Various types of postoperative care can be practiced.



* Various drainage tubes, catheters, or adhesive tapes are not included.

Accessories



Suture seal

Stoma


Electrodes

Sham tracheal cannula

Male/Female external genitalia

Practice


Various postoperative care situations are simulated.



● Tracheotomy

● Drain management
* Two holes for abdominal and thoracic drainage are prepared.

● Continuous bladder catheterization



Exchangeable parts

- Suture seal
- Stoma
- Electrodes (6 pieces/set)