

M107-3 / M107-4

Renewal

SAKAMOTO Baby (boy) SAKAMOTO Baby (girl)

- Accessories/Umbilical cord, Baby powder, Baby clothes, Simulated loose excrement (yellow), Container for mixing simulated loose excrement, Muddler for mixing simulated loose excrement, Tray



Umbilical cord
(fit-on & off available)

The facial expression became cute.



- The limbs have seamless structures.
- The head moves in every direction. (Please hold the model by supporting the head as you would do in holding a neonate.)
- The weight of the head has the same proportion as that of a neonate.
- The anterior fontanelle and ribs can be palpated.
- Water infused into the mouth flows out from the urethra so that changing diapers can be practiced.
- Simulated soft stools can be used to carry out cleaning training.

M107-5 / M107-6

Blue eyes baby has joined the lineup.



Specifications

| | |
|----------|-------------|
| Height | 49cm |
| Weight | About 3.3kg |
| Material | Silicone |

Practice

1. Bathing



2. Baby exercise and feeding



3. Dressing / undressing and changing diapers



4. Anal thermometry and various measurements



Others Usable also for the practice of umbilical treatment, cleaning of the acoustic meatuses and nostrils, and guidance at infant screenings.

Nursing

Injection

Emergency

Experience

Acupuncture