

Sakamoto Obstetric Assistant Trainer

- Weight/About 13kg ● Case size/W74×D37×H55cm
- Accessories/External genitalia (skin color/transparent), Fetus (Delivery), Lubricant oil, Baby powder

Transparent external genitalia make it possible to discern the status of fetal head.



Feature

1 Two types of external genitalia (skin color/transparent)

Can see the status of fetal head through the transparent genitalia.



2 Two types of posture

Supine position, all fours possible.

1 Obstetric assist

- Protection of the perineum



- Crowning



- Delivery of anterior shoulder



2 Can check placenta/amnion.



3 Umbilical ligation and umbilical cord cutting

Forceps can now be used as well.



▲Hemostatis using forceps



▲Umbilical cord cutting

4 Two types of posture

Supine position, all fours possible



▲All fours (Knee-chest position)

Exchangeable parts

- Exchangeable genitalia (skin color)
- Exchangeable genitalia (transparent)
- Exchangeable umbilical cord (10 pieces/set)